

IHC - Immunohistochemistry – from paraffin-embedded slides

- 1) Xylene washes (3 containers, #1, #2, #3, 3x5min) to de-paraffinize (xylene waste disposal in brown bottles under fume hood)
- 2) EtOH washes (100% x2, 95%, 70%, 50%) – 700ml 100% EtOH + 300ml H₂O from big tank (in 50ml glass slide holders), 5x5min
- 3) PBS washes, 2x 5min
- 4) Antigen retrieval/unmasking
 - a. Make buffer: 4.41g sodium citrate (our shelf) in 2l plastic container + 1l ddH₂O from big tank, read ph (rinse first with ddH₂O), bring PH to 6 by adding drops of 12N HCl using a transfer pipette, mix on shaker with a fish, add 500ml more ddH₂O
 - b. Pressure cooker: heat buffer for 2min (set time while blinking) at low pressure
 - c. Put slides in black holder, close lid, put in buffer, heat for 5min at high pressure
 - d. Release pressure, put pressure cooker bucket in sink with ice and let cool down for 15min.
- 5) Wash 2x 5min with PBS
- 6) 10% H₂O₂ (5ml hydrogen peroxide (Sheffield shelf) + 45ml PBS) – let slides sit for 20min – to inactivate endogenous peroxidase
- 7) Flick off H₂O₂, wipe slide on bottom, circle sample on top, add T-PBS (50microl Tween (improves permeability of cell membrane so that ab can enter the nucleus) in 50ml PBS = 0.1% T-PBS) – wash 3x 2min
DON'T LET SLIDES DRY OUT
- 8) Add 5% Natural Serum (against species of 2.ab) to block ??? (950microl PBS-T + 50microl serum (small brown bottle in fridge) – let slides sit for 45min.
- 9) Flick off and add avidin blocking solution – let sit for 15min.
- 10) 1x5min PBS

- 11) Flick off NGS, add a few drops **biotin blocking solution** – let it sit for **15min.**
- 12) **2x5min** PBS wash
- 13) Flick off and add **1.ab** / NO AB for NEG CONTROL: mouse active b-catenin (in fridge, red lid), 1:300, (1000microl PBS + 3microl antibody) – incubate in fridge inside slide box **overnight**
-
- 14) Warm slides for **1h** at 37C in small oven
- 15) **2x 5min** PBS wash
- 16) **2.ab**: biotinylated goat anti-mouse, 1:500, 2microl ab + 1000microl PBS-T – let sit **45min**
- 17) **2x5min** PBS wash
- 18) **ABC kit**: (avidin + biotinylated horseradish peroxidase H (HRP) – detects biotinylated 2.ab) – put a few drops of ABC reagent (grey bottle) on slides – let sit for **30min**
- 19) **2x 5min** PBS wash
- 20) **DAB solution** (peroxidase substrate to localize peroxidase) – in 15ml tube mix 5ml of ddH₂O + 2 drops buffer + 4 drops DAB (diaminobenzidine – produces brown reaction product in presence of HRP) + 2 drops HP solution
- 21) Go over to microscope, turn on in the back, wheel on 0, yellow magnification
- 22) PC: camera control: live
- 23) Add DAB solution to slides and look under the microscope for brown stain
- 24) After **5-10min** dump slides in glass jar with ddH₂O to stop the reaction
- 25) Counter-stain with **Hemotoxylin** – transfer slides for 30sec into glass jar with H (from under the hood) and then back into jar with ddH₂O, recycle H back into big glass jar, wash slides inside the glass jar with tap water to bring out blue stain (nucleus), add a little bit of ddH₂O from glass jar to slides and look under microscope
- 26) EtOH dilutions for de-hydration: add 50% EtOH into glass jar, let sit **5x5min.**, then 70%, 90%, 100%

27)Xylene 1x10min in jar#3

28)Mounting